

Conjugating CODEX® tags on antibodies of choice

1. Prepare Conjugation Reagents and Antibodies

- Purify antibodies if needed (see User Manual)
- Measure and calculate concentration of Antibodies using pre-set IgG settings of Nanodrop
- Calculate the volume of solution corresponding to **50 µg of antibody**.
- Retrieve the following reagents now:
 - Reduction Solutions 1 & 2
 - Filter Blocking Solution
- Retrieve the following reagents in ~1 hour:
 - Conjugation Solution
 - Barcodes
- Retrieve the following reagents in ~3 hours:
 - Purification Solution
 - Antibody Storage Solution

Antibody stocks:


- Purified Antibodies in PBS or in a similar buffer. Antibodies should be free of carrier proteins, such as BSA, and other chemicals such as gelatin or glycerol.

Akoya Materials

- CODEX® Conjugation Kit
 - Reduction Solution 1
 - Reduction Solution 2
 - Filter Blocking Solution
 - Conjugation Solution
 - Purification Solution
 - Antibody Storage Solution
- CODEX® Barcodes

Additional Materials

- 50kDa MWCO filter
- 1.5 ml sterile tube(s)
- 1.5 ml low bind tubes
- ddH₂O
- 0.2 ml PCR tubes
- Ice bucket
- Centrifuge for 1.5 ml tubes
- NanoDrop



2. Reduce Purified Antibody

- Label a 50kDa MWCO filter for each antibody.
- Add **500 µl** of **Filter Blocking Solution** to the top of **each 50kDa MWCO filter**.
- Spin down at **12,000g** for **2 mins**.
- Remove all liquid; discard both the liquid left on the top of the column and the flow-through solution. Use a pipette if necessary to remove all liquid on top of the column.
- Add **50 µg** of the **antibody** in a volume **100 µl** or greater.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.
- Create the **Antibody Reduction Master Mix** based on the number of CODEX® antibody conjugates.

Antibody Reduction Master Mix								
Number of Antibodies	1	2	3	4	5	6	7	8
Reduction Solution 1 [µL]	6.6	13.2	19.8	26.4	33	39.6	46.2	52.8
Reduction Solution 2 [µL]	275	550	825	1100	1375	1650	1925	2200
Total [µL]	281.6	563.2	844.8	1126.4	1378	1689.6	1971.2	2252.8

- Add **260 µl** of the **Antibody Reduction Master Mix** to the top of each filter unit. Vortex for 3 seconds.
- Incubate for **30 mins**. Exceeding 30 mins will result in irreparable damage to antibodies.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.
- Add **450 µl** of **Conjugation Solution** to the top of each column.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.

3. Prepare CODEX® Barcodes and Conjugate Antibody

- Resuspend each **Barcode** in **10 µl** of Molecular Biology Grade **Water** by pipetting up and down.
- Add **210 µl** of **Conjugation Solution** to each **suspended Barcode**. Mix by gentle pipetting. Set aside.
- Add the **CODEX® Barcode Solution** to the top of each filter. Close the lid and vortex for 3 seconds.
- Incubate for **2 hours** at **RT**

4. Purify and Collect Tagged Antibodies

- Set aside **2-5 μl** of the purified solution for QC and troubleshooting.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.
- Add **450 μl** of **Purification Solution** to the top of each column.
- Spin down tubes at **12,000 g** for **8 mins**. Discard the flow-through.
- Repeat steps c and d** two more times for **a total of three** purifications.
- After the **third centrifugation, discard** the flow-through solution.
- The top of the column will contain the remaining purified solution.
- For each antibody, **label a new tube** and lid that can hold filter units..
- Add **100 μL** of **Antibody Storage Solution** to each filter unit column.
- After it is labeled, **place** the new empty **tube upside-down** on top of the filter unit column.
- Invert the filter** for collection into the new collection tube.
- Spin solution down at **3,000g** for **2 mins**. The **final volume** in the tube should be **around 120 μl** .
- Transfer the solution to a sterile, screw-top tube for storage at **4°C** for up to **1 year**.
- Do not use these antibodies for tissue staining for at least 2 days; if used for staining sooner, you may experience high levels of background nuclear staining.

5. Quality Control of Conjugated Antibodies

- Dilute** conjugated antibodies and the unconjugated antibody (control) to a **final volume** of **13 μl** in Nuclease free water.
- Add **5 μl** of NuPAGE™ LDS Sample Buffer to each sample.
- Add **2 μl** case NuPAGE™ Sample Reducing Agent to each sample.
- Denature at **95°C** in a dry bath for **5 mins**.
- Prepare buffer by adding **25ml** of NuPAGE MOPS SDS Running Buffer and **475ml** of NuPAGE MOPS SDS Running Buffer.
- Prepare gel according to manufactures instructions.
- Pour buffer into gel tank.
- Add **5 μl** of a pre-stained protein standard 3.5-260 kDa to the gel.
- Add **20 μl** of antibody to one well each.
- Run the gel **at 200 V** for **30-40 mins** until completion.
- Place in a container filled with ddH₂O.
- Microwave** the gel until the **first bubbles form**.
- Stain the gel with Novex SimplyBlue™ SafeStain.
- Microwave** the gel again until the first **bubbles form**.
- Place the gel in a shaker for **10 mins**.
- Wash the gel with ddH₂O and leave it on the shaker until bands are visible.

Antibody:

- 1 μg of each conjugated antibody

Additional Materials

- NuPAGE™ LDS Sample Buffer (4X)
- NuPAGE™ Sample Reducing Agent (10X)
- NuPAGE™ 4-12% Bis-Tris Protein Gels
- Novex™ Sharp Pre-Stained Protein Standard - 3.5-260 kDa
- XCell SureLock™ Mini-Cell Electrophoresis System
- NuPAGE™ MOPS SDS Running Buffer (20X)
- Novex™ SimplyBlue™ SafeStain ddH₂O
- 95°C dry bath
- Microwave
- Shaker

